



Semi-Automatic Blood Pressure Monitor

Model : CH91



- Digital Technology
- Irregular Heartbeat Detection (IHB)
- 90 memories
- One-touch operation
- Latex-free patented cone cuff



Digital Technology



Irregular Heartbeat Detection (IHB)



90 memories



One-touch operation



Latex-free patented Cone cuff



AAMI

This blood pressure monitor meets the AAMI accuracy standard requirements (AAMI/ANSI/ISO 81060-2)



Cone cuff fits arm circumference:
24 - 36cm (9.4" - 14.2")
Other size cuffs available upon request

Digital Technology

Easy to use manual inflation. This manual inflate blood pressure monitor offers accurate blood pressure and pulse measurements quickly and easily. Simply inflate the cuff to the desired pressure. The cuff will automatically deflate and show your readings on the easy to read display.

Model	Qty per carton	Carton volume
CH91	18 pcs	0.051 cbm/ctn
QTY per 20'	QTY per 40'	QTY per 40'HQ
9882 pcs	19764 pcs	22932 pcs

- Storage and Transportation Condition
-10°C~60°C(14°F~140°F)
10%~90% RH, 700~1060 hPa